



## Aggregate Washers

LA-0797-02 Small | LA-0797 Large

Forney Part No.	Description
LA-0797-02	Small Aggregate Washer
LA-0797	Large Aggregate Washer

- For 115V units, plug into a grounded wall outlet.
- For 230V units, plug the input cord into the included transformer. Do not operate directly on 230 VAC.

### POWER CONNECTION

1. Plug unit into a grounded power socket.
2. Press and release the Reset button on the GFCI.
3. Verify that the red band at the base of the reset button is not visible.
4. Press the Test button to ensure the GFCI is working and power is restored.
5. Press and release the Reset button.
6. Verify that the red band is not visible.

### WARNING!

If the GFCI fails to trip when the test button is pressed, or fails to reset, the device is defective and should be returned for replacement.

If the GFCI tests properly without any power cord attached, but trips each time the appliance cord is connected, the appliance has a ground fault and needs to be repaired or replaced.



LA-0797-02 Small | LA-0797 Large

### WATER CONNECTION

The unit is supplied with 15 feet of 0.25 inch diameter tubing and one saddle valve that requires installation in a copper water line.

1. Turn off the water supply to the water line being used.
2. Drill a 15/64 inch (0.234 inch) hole in the copper water line for the saddle valve.
3. Assemble the water valve to the saddle clamping device using pipe sealant or tape.
4. Install the rubber washer over the protrusion in the saddle body and attach the assembly to the water line using the screws provided.
5. Connect the 0.25 inch tubing to the valve, turn the water on, and check for leaks.
6. Open valve by turning handle counterclockwise to the desired water flow rate.

### ACCESSORIES & REPAIR PARTS

- Small wash drum
- Large wash drum
- Vacuum lid
- 8 inch diameter sieving pan with drain
- 12 inch diameter full-height brass pan with drain
- Clear vinyl tubing
- Hose clamp
- Flexible water tube
- Rubber wheel cover

For accessory availability, replacement parts, or current part-number support, contact Forney Customer Service.

Rev: 06/2026

# FORNEY

WWW.FORNEYONLINE.COM

(800) 367-6397

2050 Jackson's Pointe Court  
Zelienople, PA 16063 USA